

GLAZING INSTRUCTIONS with DRILLING PATTERN

LEATHERDROPS





1. Download Drilling Coordinates



- Print drilling coordinates in size 1:1. (In print window set your page scaling "none".)
- <u>Alternative A:</u> Drilling with an Edger and Automatic Drilling Machine. Drill at a 90° angle to the front surface of the lens.
- <u>Alternative B:</u> Drilling by hand with drilling pattern and tools from mounting box P0027.



2. Cutting



- Measure lenses in the focimeter. Mark the central axis with a waterproof pen. Cut the lenses using a flat edge at a 1:1 scale.
- Polish edges as desired.
- Maximum lens thickness at fixing position = 6 mm (= length of BLS)



3.1. Alternative B: Drilling by Hand Sticking of Drilling Pattern on MINUS Lenses



Alternative B:

Use enclosed demo-lens as drilling pattern. Affix the demo-lens with double-adhesvie tape <u>accurate to</u> <u>contour</u> and <u>axially parallel</u> to the optical lens. Please note direction for fixing!



3.2. Alternative B: Drilling by Hand Sticking of Drilling Pattern on PLUS Lenses



Alternative B:

Use enclosed demo-lens as drilling pattern. Affix the demo-lens with double-adhesive tape <u>accurate to</u> <u>contour</u> and <u>axially parallel</u> to the optical lens. Please note direction for fixing!



3.3. Alternative B: Drilling by Hand Drilling and Deburring



- Place the correction lens on a curved pad for drilling. Drill the holes at 90° to the lens front. Drilling diameter is Ø 1,4 mm.
- 2. Deburr the holes carefully on both sides.



4. Insert BLS – Cut to Size – Expand





1. Apply protection film to avoid damage while cutting. Insert BLS into drilled holes (use a demo lens to get best fit). Cut to size, depending on the lens-thickness. Place blade parallel to the lens.

2. Expand open end of sleeves with a conical pin. This will make it easier to press in frame parts.



5. Press Frame Parts into Sleeves



Press frame parts into sleeves by hand.

Use some counter-pressure at the top of the BLS, use a demo lens as a support.



6. Fixing of Frame Parts



NOTE: Focus some <u>counter-pressure</u> with the pliers <u>on the inside of the lens</u>. Then close pliers and gently ease components together.

Use universal mounting pliers P0023 for the <u>bridge</u>, pliers P0026 for the <u>sideparts</u> (mount red sideparts AC 329 to the pliers)

TIP: Always place moveable sideparts at the top of the BLS!



7. Disassembly of the Demo Lenses



Remove demo-lens nasally by using the transparent disassembly pliers AC 212-2.



Remove demo-lens temporally by using the green disassembly pliers AC 268.