



## **Performance Insert™ a731**

Optical clip-in, suitable for all 8 and 10 base models

### **Notes on improving the patient's visual performance with the a731 clip-in** (developed in co-operation with the Institute of Ophthalmology in Germany)

To maximize the visual performance of the clip-in the following guidelines should be observed.

Sphere correction range should be +/- 4 Diopters  
Cylinder correction range should be +/- 1 Diopters

Use a precise monocular PD  
To compensate for the frame facial wrap

|              |                                |
|--------------|--------------------------------|
| Plus lenses  | Make PD 1 mm narrower each eye |
| Minus lenses | Make PD 1 mm wider each eye    |

To maximize the sports protection of these frames and the clip-ins, we recommend polycarbonate lenses.

|        |    |                                 |
|--------|----|---------------------------------|
| Box    | mm | 39 x 30 (8,69 cm <sup>2</sup> ) |
| ED     | mm | 39                              |
| Bridge | mm | 26                              |