

Silhouette®

AUTHENTIC EYEWEAR. SINCE 1964.

SPX SIGNIA

Inclination



1. Hold lens fixation firmly with mounting pliers
2. Open temple

3. Enclose temple edge firmly and adjust inclination on titanium side part into ~~requested~~ direction

Adjust Temple Opening

Reduce



Enlarge



Temple Disassembly



- Close temple
- Press it through
- Lift it out of hinge

Greasing Hinge with Teflon



After cleaning grease inside loop of hinge with
Teflon High-Precision-Lubricant P 0030.
Afterwards mount temple.

Temple Assembly



- Position titanium side part at temple
- press it against safe edge of table
- and press it in.

Temple Adjustment



- Curve temple in large moves.
- Metal needle inside temple enables „cold“ adjustment of temple end.
- If necessary slight heatening to max. 70°C possible.
- Never bend!